

英文取説納付時の注意事項

本書に含まれる以下の要素は、案件毎に修正が必要となる可能性があります。
納付時には対象製品の仕様を確認の上、必要に応じて修正を加えてください。

1) P.7 “Use specified receptacle for power source”中、電源仕様に関する表記。

<u>Electric capacity:</u>	DC401: 100V AC, 50/60Hz, 6A (When using the pump, heater: 14A) DC800: 100V AC, 50/60Hz, 11A (When using the pump, heater: 17A)
----------------------------------	---

2) P. 23 "Specification”中、電源仕様に関する表記。

3) P. 25～ "Wiring Diagram”の結線図中、入力電源表記。

4) P. 27 "Replacement Parts Table”中、各種交換部品の電圧・電気容量に関する表記。

以上



Freeze Dryer

Model

DC401/800

Instruction Manual

- First Edition -

- Thank you for purchasing "Freeze Dryer, DC Series" of Yamato Scientific Co., Ltd.
- To use this unit properly, read this "Instruction Manual" thoroughly before using this unit. Keep this instruction manual around this unit for referring at anytime.



WARNING!:

Carefully read and thoroughly understand the important warning items described in this manual before using this unit.

Yamato Scientific Co. LTD.

This paper has been printed on recycled paper.


◆ Cautions in Using with Safety	1
• Explanation.....	1
• Table of Illustrated Symbols	2
• Fundamental Matters of “WARNING!” and “CAUTION!”	3
◆ Before Using this unit	5
• Requirements for Installation.....	5
• Installation Method	8
◆ Description and Function of Each Part	9
• Main Unit	9
• Control Panel.....	10
◆ Operation Method	11
• Pre-freezing.....	11
• Procedure of Operation	12
• Defrost Operation	16
◆ Handling Precautions	17
◆ Maintenance Method	19
• Daily Inspection and Maintenance	19
◆ Long storage and disposal	20
• When not using this unit for long term / When disposing	20
◆ In the Event of Failure	21
◆ After Service and Warranty	22
◆ Specification	23
• Main Unit	23
• Optional Accessories.....	24
◆ Wiring Diagram	25
• DC401	25
• DC800	26
◆ Replacement Parts Table	27
◆ Reference	28
• List of Dangerous Substances	28


Explanation

MEANING OF ILLUSTRATED SYMBOLS

Illustrated Symbols

Various symbols are used in this safety manual in order to use the unit without danger of injury and damage of the unit. A list of problems caused by ignoring the warnings and improper handling is divided as shown below. Be sure that you understand the warnings and cautions in this manual before operating the unit.

 **WARNING!** If the warning is ignored, there is the danger of a problem that may cause a serious accident or even fatality.

 **CAUTION!** If the caution is ignored, there is the danger of a problem that may cause injury/damage to property or the unit itself.

Meaning of Symbols



This symbol indicates items that urge the warning (including the caution).
A detailed warning message is shown adjacent to the symbol.



This symbol indicates items that are strictly prohibited.
A detailed message is shown adjacent to the symbol with specific actions not to perform.



This symbol indicates items that should be always performed.
A detailed message with instructions is shown adjacent to the symbol.

Table of Illustrated Symbols

Warning



Warning,
generally



Warning,
high voltage



Warning,
high temperature



Warning,
drive train



Warning,
explosive

Caution



Caution,
generally



Caution,
electrical shock



Caution,
scald



Caution,
no road heating



Caution,
not to drench



Caution,
water only



Caution,
deadly poison

Prohibit



Prohibit,
generally



Prohibit,
inflammable



Prohibit,
to disassemble



Prohibit,
to touch

Compulsion



Compulsion,
generally



Compulsion,
connect to the
grounding
terminal



Compulsion,
install on a flat
surface



Compulsion,
disconnect the
power plug



Compulsion,
periodical
inspection

Fundamental Matters of “WARNING!” and “CAUTION!”

WARNING!



Do not use this unit in an area where there is flammable or explosive gas

Never use this unit in an area where there is flammable or explosive gas. This unit is not explosion-proof. An arc may be generated when the power switch is turned on or off, and fire/explosion may result. (Refer to page28 “List of Dangerous Substances”.)



Always ground this unit

Always ground this unit on the power equipment side in order to avoid electrical shock due to a power surge.



If a problem occurs

If smoke or strange odor should come out of this unit for some reason, turn off the power key right away, and then turn off the circuit breaker and the main power. Immediately contact a service technician for inspection. If this procedure is not followed, fire or electrical shock may result. Never perform repair work yourself, since it is dangerous and not recommended.



Do not use the power cord if it is bundled or tangled

Do not use the power cord if it is bundled or tangled. If it is used in this manner, it can overheat and fire may be caused.



Do not process, bend, wring, or stretch the power cord forcibly

Do not process, bend, wring, or stretch the power cord forcibly. Fire or electrical shock may result.



Pay special attention to the measure for flammability and handling of flammable solvent

Leaving at the temperature higher than the room temperature may vaporize the flammable material (ethanol, etc.). There might be the case that some flammable liquid might be vaporized at the temperature lower than the room temperature. The result of such careless handling could cause the fire or explosion. Do provide the vaporization with enough during the operation.



Do not disassemble or modify this unit

Do not disassemble or modify this unit. Fire or electrical shock or failure may be caused.

Fundamental Matters of “WARNING!” and “CAUTION!”

CAUTION!



During a thunder storm

During a thunderstorm, turn off the power key immediately, then turn off the circuit breaker and the main power. If this procedure is not followed, fire or electrical shock may be caused.



Do not touch ice in the trap bath

If touch the ice in the trap bath during operation or just after operation, the frostbite might be resulted from low temperature. Therefore, do not touch the ice in the trap bath.



Do not touch the cooling fin with bare hands

Do not touch the cooling fin with bare hands during maintenance, for the edge of the cooling fin is too sharp to cut your hand.

Requirements for Installation

WARNING!

1. Always ground this unit



- Connect the power plug to a receptacle with grounding connectors.
- Do not forget to ground this unit, to protect you and the unit from electrical shock in case of power surge. Choose a receptacle with grounding connectors as often as possible.
- Do not connect the grounding wire to a gas pipe, or by means of a lightning rod or telephone line. A fire or electrical shock will occur.
- Do not use a branching receptacle, which may cause the heat generation.




2. Choose a proper place for installation



- Do not install this unit in a place where:
 - ◆ Rough or dirty surface.
 - ◆ Flammable gas or corrosive gas is generated.
 - ◆ Ambient temperature below 5°C or above 30°C.
 - ◆ Ambient temperature fluctuates violently.
 - ◆ There is direct sunlight.
 - ◆ There is excessive humidity and dust.
 - ◆ There is a constant vibration.
 - ◆ Winds from the air conditioner, etc. hit the sample container directly.




- Install this unit on a stable place with the space as shown below.



3. Do not use this unit in an area where there is flammable or explosive gas



- Never use this unit in an area where there is flammable or explosive gas. This unit is not explosion-proof. An arc may be generated when the power switch is turned ON or OFF, and fire/explosion may result.
- To know about flammable or explosive gas refer to page28 "List of Dangerous Substances".




Requirements for Installation

4. Do not modify



- Modification of this unit is strictly prohibited. This could cause a failure.


Modification



5. Do not topple or tilt this unit



- Set this unit to the flattest place. This unit incorporates the refrigerator. Do not topple or tilt it. Setting this unit on rough or slope place could cause the vibration or noise, or cause the unexpected trouble or malfunction.



Requirements for Installation

CAUTION!

6. Use specified receptacle for power source



- Choose a correct power distribution board or receptacle that meets the unit's rated electric capacity.

Electric capacity: DC401: 100V AC, 50/60Hz, 6A (When using the pump, heater: 14A)
DC800: 100V AC, 50/60Hz, 11A (When using the pump, heater: 17A)

NOTE)

- ❖ There could be the case that the unit does not run even after turning ON the power. Inspect whether the voltage of the main power is lowered than the specified value, or whether other device(s) uses the same power line of this unit. If the phenomena might be found, change the power line of this unit to the other power line.
- ❖ Starburst connection with a branching receptacle or extended wiring with a cord reel lowers electrical power voltage, which may cause the degradation of refrigeration capability.

7. Before/after installing



- It may cause injury to a person if this unit falls down or moves by the earthquake and the impact. etc..To prevent, take measures that the unit cannot fall down, and not install to busy place.
- Though this unit has the air-cooled refrigerator, the device exhausts the heat. Do provide the vaporization with enough so as not to raise the ambient temperature caused by the exhaust of the heat, or install this unit with its air controlled completely. If the ambient temperature becomes high, the operation efficiency becomes worse, and could cause the malfunction of the device by high temperature and humidity.

8. Handling of power code





- Do not entangle the power cord. This will cause overheating and possibly a fire.
- Do not bend or twist the power cord, or apply excessive tension to it. This may cause a fire and electrical shock.
- Do not lay the power cord under a desk or chair, and do not allow it to be pinched in order to prevent it from being damaged and to avoid a fire or electrical shock.
- Keep the power cord away from any heating equipment such as a room heater. The cord's insulation may melt and cause a fire or electrical shock.



- If the power cord becomes damaged (wiring exposed, breakage, etc.), immediately turn off the power at the rear of this unit and shut off the main supply power. Then contact your nearest dealer for replacement of the power cord. Leaving it may cause a fire or electrical shock.
- Connect the power plug to the receptacle which is supplied appropriate power and voltage.

Installation Method

1. Unlock the stopper of the caster.
Pulling up the lever of the stopper for caster releases the lock.
(Only the two casters in front of the unit are attached the stopped.)
2. Move the device to the place to be installed.
- ❖ If there is a step on the floor, the too strong impact is given to the caster, and could give the damage. In that case, move the device by lifting at the step.
3. When the installation place is determined, pull down the lever of the stopper for caster, and lock them.
4. Drain Tap Check
 - Check whether the drain tap is detached or not.




5. Power Plug Connection
Check the power of the earth leakage breaker is turned "OFF", and plug the power cord in the receptacle.
6. Prepare the vacuum pump that the exhaust amount is 50ℓ /min or more, and the ultimate pressure is 1.0×10^{-1} Pa. This pump must also have a check-valve.
7. Install the vacuum valve to all of the valve mounting holes on the chamber.
8. Prepare a spare freeze dryer that can retain temperature sufficiently lower than the eutectic (crystal) point of the sample.


Description and Function of Each Part

Main Unit

Front view




Rear view



Description and Function of Each Part

Control Panel



Part Name		Function
①	Refrigerator Switch :	The refrigerator activates by pressing "ON" of this switch.
②	Vacuum Pump Switch :	Enables pump operative.
③	Pump Start Switch :	Starts pump operation.
④	Pirani Gauge :	Displays degree of vacuum.
⑤	Defrost Switch :	Melts frozen ice in the trap bath.
⑥	Refrigerator Trouble Lamp :	This lamp lights on when the refrigerator is in overload.

Pre-freezing

At pre-freezing process, quickly and completely freeze to the inside at sufficient lower temperature than the eutectic (crystal) point* of the sample. Make the ice layer as thin as possible because thick ice layer might be easily melted by drying. Moreover, keep in mind that it becomes easy to receive the heat of external temperature, and it may cause melting-ice if it freezes widely on the container surface when the eutectic point of the sample is low.

Keep in mind that if there are too many throughput, or if the eutectic point of the sample is low, a melting-ice phenomenon arises and this may cause bumping or scattering of the sample.

When installing the container and starting freeze-drying after pre-freezing process, set 3-4 minutes' interval till next sample drying starts. (This differs depending on the container and sample types.)

Melting-ice may be caused if the interval is too short.

Eutectic (Crystal) Point



- Although pure water is frozen at 0 Celsius degree, solution is not frozen immediately below 0 Celsius degree, and the solution becomes soft ice condition because the frozen pure water and the concentration solution that is not frozen are mixed.
Furthermore, if freezing is continued, all the solution will freeze uniformly under containing its solutes, and will become solid ice in whole.
This temperature is called as the eutectic (crystal) point of the solution.

Procedure of Operation





Before operation, prepare a drying chamber, a vacuum pump, a vacuum hose, a tandem tube, etc.

Vacuum Pump: PD52 (YAMATO Product)





Exhaust amount: 50ℓ /min, Ultimate pressure: 6.7×10^{-2} Pa,
Diameter of inlet pipe: 18mm (or equivalent)

<p>1. Supply the power</p> 	<ul style="list-style-type: none"> • Turn "ON" the earth leakage breaker.
<p>2. Operation of the refrigerator</p> <div style="display: flex; justify-content: space-around;"> <div data-bbox="193 779 384 999"> <p>DEFROST</p> <p>ON</p> <p>OFF</p> </div> <div data-bbox="416 779 600 999"> <p>REFRIGE</p> <p>ON</p> <p>OFF</p> </div> </div>	<ul style="list-style-type: none"> • Turn "ON" the refrigerator switch. The refrigerator activates. ❖ There is the case that the starting sound of the refrigerator might become higher depending on the operation initial status and ambient temperature. However, this phenomenon is not abnormal.
<p>3. Installation of chamber</p> 	<ol style="list-style-type: none"> 1. Make sure that there is no dust adhered to the trap packing or the chamber for flask installation. <ul style="list-style-type: none"> ❖ If any dent or dust is existed, degree of vacuum in the trap bath becomes low. 2. Gently put the chamber for flask installation etc. on trap packing. <ul style="list-style-type: none"> ❖ If it has not stuck to trap packing enough, it will become the cause of a vacuum leak. ❖ The same is said of the case of a manifold etc.




Procedure of Operation

<p>4. Connection with the vacuum pump</p> 	<ul style="list-style-type: none"> • Connect the vacuum hose to the exhaust port on the left side. (Exhaust port: outer shape 17 mm) (Inside dia. of applicable vacuum hose: 12 X outer shape 30 mm) Connect the opposite side of the vacuum hose to the inlet of the vacuum pump. ❖ The vacuum hose is not included. ❖ Seal strength shall be improved by applying the vacuum seal (silicon grease compound manufactured by TORAY, H.V.G, etc) inside of the vacuum hose at connecting.
<p>5. Connect the vacuum pump with the rear receptacle</p> 	<ul style="list-style-type: none"> • Insert the power code of the vacuum pump into the 6A receptacle of the rear side. • Keep the switch of the vacuum pump ON. ❖ Use the receptacle below 8A in total.
<p>6. Make all valves for "VENT" to turn up</p> 	<ul style="list-style-type: none"> • Make all valves for "VENT" face to turn up, and close the path of the chamber side.
<p>7. Turn ON the vacuum pump switch. Then, press Start switch once</p> 	<ul style="list-style-type: none"> • After the trap bath becomes cold enough (if the lowest temperature attainment time is over), start the vacuum pump operation. • The Lowest temperature attainment time: (Outside temperature: 20°C, no load) DC401: 50min., DC800: 80min. This is guideline. (Condensation might be caused at the bottom of the chamber.)





Procedure of Operation

<p>8. Check the vacuum gauge</p> 	<ul style="list-style-type: none"> Make sure that inside of the chamber attains vacuum with the vacuum gauge. (10Pa-20Pa)
<p>9. Install a sample container, then, turn up "VAC" face of the vacuum valve</p> 	<ul style="list-style-type: none"> Install the pre-frozen sample container to the vacuum valve, then, turn up "VAC" face of the vacuum valve, and open the path of the chamber side. The sample container is decompressed and drying process starts.
<p>10. Container installation after second piece</p> 	<ul style="list-style-type: none"> When switching the vacuum valve setting from "VENT" to "VAC", inside of the chamber is pressurized. In case of multiple samples processing at one time, wait (approx. 3-4 min.) till the pressure is returned, and then, switch the next vacuum valve from "VENT" to "VAC". ❖ If installing multiple samples at one time, the degree of vacuum decreases and the pre-frozen sample is melted. This may cause the sample container broken.
<p>11. End of freeze-drying</p> 	<ul style="list-style-type: none"> Finish the freeze-drying process after checking the sample condition.

Procedure of Operation


<p>12. After freeze-drying process finished, turn "VENT" face of the vacuum valve upward</p> 	<ul style="list-style-type: none"> • Turn "VENT" face of the vacuum valve upward in order to return the inside of the sample container to normal atmosphere pressure. Then, remove the container. Keep the container securely by hands so as not to fall.
<p>13. Return inside of the chamber to normal atmosphere pressure</p> 	<ul style="list-style-type: none"> • When drying of all samples are completed, and the container is removed, turn one of the vacuum valves gradually so that "VAC" face turns upward, and return the inside of the chamber to normal atmosphere pressure.
<p>14. Turn OFF the vacuum pump</p> 	<ul style="list-style-type: none"> • Switch OFF the vacuum pump switch after the inside of the chamber is returned to normal atmosphere pressure.

Defrost Operation


<p>1. Remove the chamber for flask installation</p>  <p>Chamber for flask installation (Optional)</p>	<ul style="list-style-type: none"> Remove the chamber for flask installation. If ice coating is formed on trap bath wall, turn ON the defrost switch.
<p>2. Remove ice</p> 	<ul style="list-style-type: none"> If the ice of the ice coating side is melted by defrosting, take out the ice and turn OFF the defrost switch. Defrosting is auto-stopped. DC401: Approx. 20min. DC800: Approx. 45min. ❖ Switch OFF the defrost switch by manual.
<p>3. Drain water</p>  <p>Drain tap</p>	<ul style="list-style-type: none"> Remove the drain tap, and drain the water remained in the bath with the drainage hose.
<p>4. Stop operation</p> 	<ol style="list-style-type: none"> 1. Turn "OFF" the refrigerator switch. 2. Turn "OFF" the earth leakage breaker.

WARNING!

If a problem occurs


-  If smoke or strange odor should come out of this unit for some reason, turn off the power key right away, and then turn off the breaker and the main power. Immediately contact a service technician for inspection. If this procedure is not followed, fire or electrical shock may result. Never perform repair work yourself, since it is dangerous and not recommended.

Measure for flammability and handling of flammable solvent


-  This unit is not designed as the explosion-proof construction. Pay special attention to the handling of the sample to be handled with this unit on the consumption with the explosive material, flammable material, and similar ones. The flammable material may be vaporized by leaving it at the temperature higher than room temperature, and could cause the fire or explosion. When handling such material, provide ventilation with enough before the operation. (Refer to page 28 "List of Dangerous Substances".)

CAUTION!


Maximum trap capacity

-  The max. trap capacity of the bath in this DC401/800 is approx. 0.6ℓ /1.0ℓ . Excessive amount shall cause significant unit performance degradation.


Do not step on this unit

-  Do not step on this unit. It will cause injury if this unit fall down or break.


Do not put anything on this unit

-  Do not put anything on this unit. It will cause injury if fall.


During a thunder storm

-  During a thunderstorm, turn off the power key immediately, then turn off the circuit breaker and the main power. If this procedure is not followed, fire or electrical shock may be caused.


Countermeasure for stop operation during night or long-term stop

-  In case of stopping operation during night or long-term, toggle the power switch to "OFF".


Power supply voltage

-  The power supply voltage must be within +/-10% of the rated voltage.

Do not use solutions

-  The rubber packing of the unit and U-ring on the chamber cap for the vacuum valve and the flask are made of chloroprene rubber. Keep in mind that they may be immersed by acid, halogen, aromatic series, ester, and the oxo solvent.

Restarting

-  Stop operation and wait over 5 minutes, and then, restart operation.

Upper limit of inside bath temperature



Do not operate continuously when the temperature of the bath is higher than -20°C.

Vacuum grease application



If there is any dirt adhered to each connecting part of the vacuum system, the degree of vacuum might be decreased. In this case, wipe the dirt off, and then, apply the vacuum grease.

In case that the stem tube of the vacuum valve is hard to rotate, remove the stem tube once and apply the vacuum grease, and then install the stem tube again.

When returning to normal atmosphere pressure



Before turning OFF the vacuum pump, make sure that inside of the chamber is returned to normal atmosphere pressure after operation finished. If turning OFF the vacuum pump under vacuum condition, the oil inside of the pump might be flown back in the chamber. (Refer to “12.

After freeze-drying process finished, turn “VENT” face of the vacuum valve upward” in P.15 for the procedure of vacuum relief.)

Dimension of ampoule



Use an ampoule that diameter of the connecting hole is approx. 7-9mm since the inner diameter of the ampoule adapter is 7mm.

Eggplant flask



Use TS29 for grinding eggplant flask.


Daily Inspection and Maintenance

WARNING!

- Disconnect the power cable from the power source when doing an inspection or maintenance unless needed.
- Perform the daily inspection and maintenance after returning the temperature of this unit to the normal one.
- Do not disassemble this unit.
- Do not touch the cooling fin with bare hands.


CAUTION!

- Use a well-drained soft cloth to wipe dirt on this unit. Do not use benzene, thinner or cleanser for wiping. Do not scrub this unit. Deformation, deterioration or color change may result in.




Monthly maintenance


- Check the earth leakage breaker function.
 1. Connect the power cord.
 2. Turn the breaker on.
 3. Push the red test switch by a ballpoint pen etc.
 4. If there is no problem, the earth leakage breaker will be turned off.



Cleaning of cooling fin

- Clogging of the cooling fin could cause the deterioration of the cooling performance, and also cause the malfunction of the refrigerator. The clogged status differs depending on the surrounding condition or operation time. Clean the cooling fin periodically.
- ❖ Loosen the mounting screws (4 screws) of the ventilation port cover, remove the cover of the ventilation port, and remove the dust attached to the surface of the cover using the vacuum cleaner.
- ❖ After cleaning the cooling fin, attach it in inverse procedure.

 Take care not to crush the fin during cleaning.



For any questions, contact the dealer who you purchased this unit from, or the nearest sales division in our company.

Long storage and disposal

When not using this unit for long term / When disposing

CAUTION!

When not using this unit for long term...

- Turn off the power and disconnect the power cord.

WARNING!

When disposing...

- Keep out of reach of children.

Environmental protection should be considered

We request you to disassemble this unit as possible and recycle the reusable parts considering to the environmental protection. The feature components of this unit and materials used are listed below.

Component Name	Material			
Parts of Main Unit				
Casing	Bonderizing steel plate baked with melamine resin coating			
Inner bath, Cover	Stainless steel SUS304, Acryl			
Production plates	Polyester (PET) resin film			
Corner	Alkylbenzenesulfied (ABS) resin			
Trap packing	Silicon rubber			
Electrical Parts				
Switches, Relays	Composites with resin and others			
Power cord & wiring materials and others	Composites with synthetic rubber, copper, nickel and others			
Parts of Refrigeration System				
Refrigerator unit	Iron, Synthetic oil and others			
Cooling fin	Aluminum, Copper			
Parts of Piping				
Drain hose	Natural rubber			
Drain tap	Polyacetal resin			
Drain tap holder	Bonderizing steel plate baked with melamine resin coating			
Pipe cover	Polyurethane sponge			
Pipe	Copper			
Sealed Cooling Medium for Refrigerator				
Cooling medium		DC401	DC800	Ask the specialist for the dealing of cooling medium.
	R404A	225g	450g	
	R23	-	53ℓ	

Trouble Shooting

Condition	Possible Causes
The device does not start when turning on the power switch.	<ul style="list-style-type: none">• Power plug is not connected to the receptacle correctly.• Power failure.• Earth leakage breaker is turned to "OFF"
Not fallen the temperature.	<ul style="list-style-type: none">• The cooling fin is clogged.• The cooling sample is overheated.• The ambient temperature is exceeding 30°C.• The peripheral of the ventilating port is shut down.

In lighting on the lamp,

Error Sign	Cause	Remedy
Turned on the refrigerator trouble lamp	Overload of refrigerator	<ul style="list-style-type: none">• Immediately turn off the power, remove the cause of the error referring to the "Trouble Shooting" (Not fallen the temperature) above, and turn on the power again after passing the certain time.• In case of lighting the alarm lamp again, make a call to the service office.

In the case if the error other than listed above occurred, turn off the power switch and primary power source immediately. Contact the shop of your purchase or nearest Yamato Scientific Service Office.

In Case of Request for Repair

If the failure occurs, stop the operation, turn OFF the power switch, and unplug the power plug. Please contact the sales agency that this unit was purchased, or the Yamato Scientific's sales office.

< Check following items before contact >

- ◆ Model Name of Product
 - ◆ Production Number
 - ◆ Purchase Date
 - ◆ About Trouble (in detail as possible)
- } See the production plate attached to this unit.

Minimum Retention Period of Performance Parts for Repair

The minimum retention period of performance parts for repair of this unit is 7 years after discontinuance of this unit.

The "performance part for repair" is the part that is required to maintain this unit.

Main Unit

Model		DC401	DC800
Method		Storage water heater body cooling method	
Performance	Dehumidifying capacity	Max. 0.6kg	Max. 1.0kg
	Lowest reached temperature ※	-45°C	-80°C
	Time required to reach lowest temperature ※	50min. (20°C to -45°C)	80min. (20°C to -80°C)
Configuration	Refrigerator	Air cooling, 400W	Air cooling, 350W × 2
	Cooling medium	R404a	R404a/R23
	Pressure indication gauge	Pirani gauge	
	Exhaust port	φ 17mm	
	Chamber material	SUS304, Cylindrical	
	Dimensions of chamber	φ 153 × H235mm	
	Capacity	Approx. 4L	
	Usable ambient temperature	5 to 30°C	
Standards	External dimensions	W340 × D450 × H720 mm	W405 × D455 × H850 mm
	Weight	Approx. 43Kg	Approx. 54Kg
	Power supply	100V AC, 50/60Hz, 6A	100V AC, 50/60Hz, 11A
Accessories		Vacuum silicon grease (TORAY H.V.G), Vacuum hose (φ 12 × φ 30 × 1.5m), Instruction manual for main unit, Instruction manual for Pirani gauge	


※: Operation conditions: Outside air temperature = 20°C, no load

Optional Accessories

Name	Specification		Product code
Chamber for flask installation	Mounting cock	Inner dia.: 18.5mm	212560
	Mounting cock pitch	96mm	
	Number of port	12	
	Dimensions	φ 195 × H303 mm	
Manifold A	Mounting cock	Inner dia.: 18.5mm	212561
	Mounting cock pitch	80mm	
	Number of port	8	
	Dimensions	W304 × D60 × H263 mm	
Manifold B	Mounting cock	Inner dia.: 18.5mm	212562
	Mounting cock pitch	80mm	
	Number of port	16	
	Dimensions	W624 × D60 × H263 mm	
Dry chamber	Number of shelves	1	212563
	60mm dish, Number of falls	7	
	Temperature adjustment	30°C ± 2°C ※	
	Dimensions	φ 252 × H240 mm	
Stopper type dry chamber	Number of shelves	1	212564
	60mm dish, Number of falls	7	
	Temperature adjustment	30°C ± 2°C ※	
	Dimensions	φ 252 × H425 mm	
Dry flask	120mℓ ・ 5 pcs.		212820
	250mℓ ・ 5 pcs.		212821
Dry flask cap (with glass joint)	Straight ・ 5 pcs.		212570
	45° bent ・ 5 pcs.		212571
Serum bottle (Vial bottle)	50mℓ ・ 10 pcs.		212814
	30mℓ ・ 10 pcs.		212815
	10mℓ ・ 10 pcs.		212816
Adapter for ampoule	Single ・ 5 pcs.		212572
	Double ・ 5 pcs.		212573
	Triple ・ 5 pcs.		212574
Adapter for test tube (with glass joint)	Straight		212590
	45° bent		212591
Adapter for egg plant flask (female)	TS29 grinding or equivalent		212569
Adapter for egg plant flask (male)	TS29 grinding		212597
Glass joint	Straight		212598
	45° bent		212599
Micro tube holder	1.5mℓ × 24 pcs. mountable		212580


※: Conditions: Room temperature = 20°C, no load

DC401



ELB	Earth leakage breaker	X3	Relay (Refrigerator)	PAsw	Switch (Vacuum gauge)	CN	Receptacle
P	Terminal block	X4	Relay (Refrigerator stop)	PSsw	Pressure switch	T1	Timer (Auto defrost)
FM	Fan motor (Refrigerator)	X5	Relay (Pressure error)	VAC	Lamp (VACsw built-in)	T2	Timer
MV1	Solenoid valve (Refrigerator)	X6	Relay (Solenoid valve)	DEF	Lamp (DEFsw built-in)	PTC-ASM	Start relay
MV2	Solenoid valve (Defrost)	VACsw1	Switch (Pump)	REF	Lamp (REFsw built-in)	OVR	Overload relay
MV3	Solenoid valve (VAC)	VACsw2	Switch (Pump start)	PS	Lamp (Pressure)	C1	Operation condenser
X1	Relay (VAC)	DEFsw	Switch (Defrost)	PA	Pirani gauge	C2	Start condenser
X2	Relay (Refrigerator, Pump stop)	REFsw	Switch (Refrigerator)	TH	Sensor (Vacuum gauge)	RF	Refrigerator

DC800



ELB	Earth leakage breaker	X4	Relay (Refrigerator stop)	VAC	Lamp (VACsw built-in)	T3	Timer (Refrigerator delay)
P	Terminal block	X5	Relay (Pressure error)	DEF	Lamp (DEFsw built-in)	PTC	Start relay
FM	Fan motor (Refrigerator)	X6	Relay (Solenoid valve)	REF	Lamp (REFsw built-in)	OVR1, 2	Overload relay
MV1	Solenoid valve (Refrigerator)	VACsw1	Switch (Pump)	PS	Lamp (Pressure)	C1, 3	Operation condenser
Mv2	Solenoid valve (Defrost)	VACsw2	Switch (Pump start)	PA	Pirani gauge	C2, 4	Start condenser
Mv3	Solenoid valve (VAC)	DEFsw	Switch (Defrost)	TH	Sensor (Vacuum gauge)	RF1, 2	Refrigerator
X1	Relay (VAC)	REFsw	Switch (Refrigerator)	CN	Receptacle		
X2	Relay (Refrigerator, Pump stop)	Pasw	Switch (Vacuum gauge)	T1	Timer (Auto defrost)		
X3	Relay (Refrigerator)	PSsw	Pressure switch	T2	Timer		

Replacement Parts Table

Common Parts

Symbol	Part Name	Specification	Manufacturer	Code No.
MV1	Solenoid valve (Refrigerator)	NEV-603DXF	Saginomiya	3-02-006-0004
MV2	Solenoid valve (Defrost)	SEV-502DXF	Saginomiya	3-02-006-0003
MV3	Solenoid valve (VAC)	AG33022 AC100V	CKD	3-20-001-0018
ELB	Earth leakage breaker	NV-L22GR 20A	Mitsubishi	LT00029776
X3, 4, 6	Relay	LY1F AC100V	OMRON	LT-00000992
X1, 2	Relay	LY3F AC100V	OMRON	LT00000993
T2, 3	Timer (Auto defrost)	ADX11134	Matsushita	2-05-000-0053
T1	Timer	ADX11184	Matsushita	LT00000994
DEFsw PEFsw	Switch	CW-SB21NMKZMEF	Nihon Kaiheiki	2-55-000-0011
VACsw2	Switch	CW-SB21NYKZYEF	Nihon Kaiheiki	2-55-000-0017
VACsw1	Switch	A3CT-90A0-Y	OMRON	LT00000995
PS	Lamp (Pressure)	BN5665L AC100V	Satoh Parts	LT00021961
CN	Receptacle	AC-R02MB12	ECHO	LT00033205
PA	Pirani gauge	GP-1 (WP-02 with sensor)	ULVAC	LT00001004
PSsw	Pressure switch	VHP-F	Fuji Kohki	3-18-000-0006

For DC401

Symbol	Part Name	Specification	Manufacturer	Code No.
RF	Compressor	RL4075HA	Hitachi	LT00028782
FM	Fan motor	SE4-C041NP	Sanyo C&C	3-01-006-0006
P	Terminal block	MKH-250ABC-4P	Terminal	LT00035672
X5	Relay	LY2F AC100V	OMRON	2-05-000-0035

For DC800

Symbol	Part Name	Specification	Manufacturer	Code No.
RF1	Compressor	C-2SN350LOU	Sanyo C&C	3-01-006-0005
RF2	Compressor	C-2SN350LOU (oil less)	Sanyo C&C	3-01-006-0008
FM	Fan motor	SE4-D11LP	Sanyo C&C	3-01-006-0014
P	Terminal block	MKH-250ABC-10P	Terminal	LT00035676
X5	Relay	LY3F AC100V	OMRON	LT00000993

List of Dangerous Substances



Never use explosive substances, flammable substances and substances that include explosive or flammable ingredients in this unit.

EXPLOSIVE

EXPLOSIVE:	Ethylene glycol dinitrate (nitro glycol), Glycerin trinitrate (nitroglycerine), Cellulose nitrate (nitrocellulose), and other explosive nitrate esters
	Trinitrobenzene, Trinitrotoluene, Trinitrophenol (picric acid), and other explosive nitro compounds
	Acetyl hidroperoxide (peracetic acid), Methyl ethyl ketone peroxide, Benzyl peroxide, and other organic peroxides

FLAMMABLE

IGNITING:	Lithium (metal), Potassium (metal), Sodium (metal), Yellow phosphorus, Phosphorus sulfide, Red phosphorus, Celluloid compounds, Calcium carbide, Lime phosphate, Magnesium (powder), Aluminum (powder), Powder of metals other than magnesium and aluminum, Sodium hydrosulfite
OXIDIZING:	Potassium chlorate, Sodium chlorate, Ammonium chlorate, and other chlorate
	Potassium perchlorate, Sodium perchlorate, Ammonium perchlorate, and other perchlorate
	Potassium peroxide, Sodium peroxide, Barium peroxide, and other inorganic peroxide
	Potassium nitrate, Sodium nitrate, Ammonium nitrate, and other nitrate
	Sodium chlorite and other chlorites
	Calcium hypochlorite and other hypochlorites
INFLAMMABLE LIQUID:	Ethyl ether, Gasoline, Acetaldehyde, Propylene chloride, Carbon disulfide, and other flammable substances having a flash point of lower than -30°C
	Normal hexane, ethylene oxide, acetone, benzene, methyl ethyl ketone, and other flammable substances having a flash point of -30°C or higher but lower than 0°C
	Methanol, Ethanol, Xylene, Pentyl acetate (amyl acetate), and other flammable substances having a flash point of 0°C or higher but lower than 30°C
	Kerosene, Light oil (gas oil), Oil of turpentine, Isopentyl alcohol (isoamyl alcohol), Acetic acid, and other flammable substances having a flash point of 30°C or higher but lower than 65°C
FLAMMABLE GAS:	Hydrogen, Acetylene, Ethylene, Methane, Propane, Butane, and other flammable substances which assume a gaseous state at 15°C and 1 atm

(Source: Appendix Table 1 of Article 6 of the Industrial Safety and Health Order in Japan)

Responsibility

Please follow the instructions in this document when using this unit. Yamato Scientific has no responsibility for the accidents or breakdown of device if it is used with a failure to comply. Never conduct what this document forbids. Unexpected accidents or breakdown may result in.

Note

- ◆ The contents of this document may be changed in future without notice.
- ◆ Any books with missing pages or disorderly binding may be replaced.

Instruction Manual for

Freeze Dryer

Model DC401/800

First Edition Aug. 7, 2009

Revision May. 30, 2012

Yamato Scientific Co., Ltd.

2-1-6 Nihonbashi Honcho, Chuo-ku,
Tokyo, 103-8432, Japan

<http://www.yamato-net.co.jp>